

Supplemental File 1. BEAPP Inputs for figures 2, 3, and 4

	Resting								Auditory
	Raw	PREP	Filter	HAPPE	CSD	REST	MskArt	RejPostSeg	ITPC
grp_proc_info.src_dir	{'/Volumes/ISP/Resting'}; ¹								{'/Volumes/ISP/Auditory'}; ¹
grp_proc_info.beapp_curr_run_tag	'Raw';	'PREP';	'Filt';	'HAPPE';	'CSD';	'REST';	'MskArt';	'RejPostSeg';	'ITPC';
grp_proc_info.beapp_prev_run_tag ²	";	Raw';	Raw';	Raw';	PREP';	'PREP';	PREP';	PREP';	";
grp_proc_info.beapp_advinputs_on	0;								0;
grp_proc_info.beapp_toggle_mods{'format',{'Module_On','Module_Export_On'}}	[1,1];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[1,1];
grp_proc_info.beapp_toggle_mods{'prepp',{'Module_On','Module_Export_On'}}	[0,0];	[1,1];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[1,1];
grp_proc_info.beapp_toggle_mods{'filt',{'Module_On','Module_Export_On'}}	[0,0];	[0,0];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];
grp_proc_info.beapp_toggle_mods{'rsamp',{'Module_On','Module_Export_On'}}	[0,0];	[0,0];	[0,0];	[1,1];	[0,0];	[0,0];	[1,1];	[1,1];	[1,1];
grp_proc_info.beapp_toggle_mods{'ica',{'Module_On','Module_Export_On'}}	[0,0];	[0,0];	[0,0];	[1,1];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];
grp_proc_info.beapp_toggle_mods{'rereference',{'Module_On','Module_Export_On'}}	[0,0];	[0,0];	[0,0];	[0,0];	[1,1];	[1,1];	[0,0];	[0,0];	[1,1];
grp_proc_info.beapp_toggle_mods{'detrend',{'Module_On','Module_Export_On'}}	[0,0];	[0,0];	[0,0];	[0,0];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];
grp_proc_info.beapp_toggle_mods{'segment',{'Module_On','Module_Export_On'}}	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];
grp_proc_info.beapp_toggle_mods{'psd',{'Module_On','Module_Export_On'}}	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[1,1];	[0,0];
grp_proc_info.beapp_toggle_mods{'itpc',{'Module_On','Module_Export_On'}}	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[0,0];	[1,1];
FORMAT									
grp_proc_info.src_format_typ	1;								2;
grp_proc_info.src_presentation_software									1;
grp_proc_info.src_data_type	1;								2;
grp_proc_info.src_linenoise	60;								60;
grp_proc_info.src_unique_nets	{ 'HydroCel GSN 128 1.0','Geodesic Sensor Net 64 2.0'};								{ 'HydroCel GSN 128 1.0', 'Geodesic Sensor Net 64 2.0'};
grp_proc_info.epoch_inds_to_process	[];								[];
grp_proc_info.src_eeg_vname	{ 'Category 1 Segment1','Category 1','Category1'};								
grp_proc_info.event_tag_offsets	0;								'input_table'; ³
grp_proc_info.behavioral_coding.events									
grp_proc_info.behavioral_coding.keys									
grp_proc_info.behavioral_coding.bad_value									
PREP									
grp_proc_info.beapp_toggle_mods{'prepp','Module Xls Out On'}		1;							1;
FILTER									
grp_proc_info.beapp_filters{'Notch','Filt On'}			1;	0;	0;	0;	0;	0;	0;
grp_proc_info.beapp_filters{'Lowpass','Filt On'}			1;	1;	1;	1;	1;	1;	1;
grp_proc_info.beapp_filters{'Lowpass','Filt Cutoff Freq'}			80;	100;	100;	100;	100;	100;	100;
grp_proc_info.beapp_filters{'Highpass','Filt On'}			1;	1;	1;	1;	1;	1;	1;
grp_proc_info.beapp_filters{'Highpass','Filt Cutoff Freq'}			4;	1;	1;	1;	1;	1;	1;
grp_proc_info.beapp_filters{'Cleanline','Filt On'}			0;	0;	0;	0;	0;	0;	0;
RESAMPLING									
grp_proc_info.beapp_rsamp_srate				250;			250;	250;	250;
INDEPENDENT COMPONENTS ANALYSIS									
grp_proc_info.beapp_ica_type				2					
grp_proc_info.beapp_toggle_mods{'ica','Module Xls Out On'}				1;					
grp_proc_info.happe_additional_chans_lbls{1}				[13,112];					
grp_proc_info.happe_additional_chans_lbls{2}				[9,58];					
REREFERENCING									
grp_proc_info.reref_typ					2;	4;			2;
grp_proc_info.beapp_reref_chan_inds									
DETRENDING									
grp_proc_info.dtrend_typ					1;	1;	1;	1;	1;

	Resting								Auditory
	Raw	PREP	Filter	HAPPE	CSD	REST	MskArt	RejPostSeg	ITPC
SEGMENTING									
grp_proc_info.beapp_toggle_mods{'segment','Module_Xls_Out_On'}									
grp_proc_info.art_thresh				40;	3000;	100;	100;	100;	3000;
grp_proc_info.beapp_reject_segs_by_amplitude	0;	0;	0;	0;	0;	0;	0;	1;	1;
grp_proc_info.beapp_happe_segment_rejection	0;	0;	0;	1;	0;	0;	0;	0;	0;
grp_proc_info.segment_linear_detrend	0;	0;	0;	0;	0;	0;	0;	0;	0;
grp_proc_info.win_select_n_trials									
grp_proc_info.beapp_baseline_msk_artifact	0;	0;	0;	0;	1;	1;	1;	0;	0;
grp_proc_info.beapp_baseline_rej_perc_above_threshold									
grp_proc_info.win_size_in_secs	1;	1;	1;	1;	1;	1;	1;	1;	
grp_proc_info.beapp_event_code_onset_strs									{'stm+'};
grp_proc_info.beapp_event_eprime_values.condition_names									{'Standard', 'Native', 'Non-Native'};
grp_proc_info.beapp_event_eprime_values.event_codes(:,1)									[1,2,3];
grp_proc_info.beapp_event_eprime_values.event_codes(:,2)									[10,12,13];
grp_proc_info.beapp_event_eprime_values.event_codes(:,3)									[11,12,13];
grp_proc_info.evt_seg_win_start									-0.1;
grp_proc_info.evt_seg_win_end									0.8;
grp_proc_info.evt_analysis_win_start									
grp_proc_info.evt_analysis_win_end									
grp_proc_info.evt_trial_baseline_removal									
grp_proc_info.evt_trial_baseline_win_start									
grp_proc_info.evt_trial_baseline_win_end									
POWER									
grp_proc_info.bw(1,1:2)									
grp_proc_info.bw_name(1)									
grp_proc_info.bw(2,1:2)									
grp_proc_info.bw_name(2)									
grp_proc_info.bw(3,1:2)									
grp_proc_info.bw_name(3)									
grp_proc_info.bw(4,1:2)									
grp_proc_info.bw_name(4)									
grp_proc_info.bw(5,1:2)									
grp_proc_info.bw_name(5)									
grp_proc_info.bw_total_freqs									
grp_proc_info.psd_win_typ	1;								
grp_proc_info.psd_interp_typ	1;								
grp_proc_info.beapp_toggle_mods{'psd','Module_Xls_Out_On'}	1;								
INTER TRIAL PHASE COHERENCE									
grp_proc_info.beapp_itpc_params.win_size									0.256;
grp_proc_info.beapp_toggle_mods{'itpc','Module_Xls_Out_On'}									1;

Grayed out cells are where inputs can be left as default values, since relevant modules are not being run. See user inputs file in BEAPP, and user guide, for further information. Users will need to set grp_proc_info.src_dir to the folder where their EEG files to be run are located. ²Files were run in the order specified in the table from left to right, so that ‘format’ and ‘prepp’ modules did not need to be repeated for subsequent runs. ³Because an offset table is specified for auditory data, offsets must be specified in beapp_file_info_table.mat.

Information table for baseline (.mat) files:

Current Folder

Name ▾

beapp

set_beapp_path.m

set_beapp_def.m

prepare_to_run_main.m

beapp_main.m

beapp_gui.m

beapp_configure_settings.m

user_inputs

rerun_fselect_table.mat

beapp_userinputs.m

beapp_set_input_file_locations.m

beapp_file_info_table.mat

beapp_advinputs.m

run_templates

reference_data

Packages

functions

documentation

Variables – beapp_file_info_table

beapp_file_info_table

10x3 table

	1	2	3	4
	FileName	SamplingRate	NetType	
1	'baselineEEG01.mat'	250	'HydroCel GSN 128 1.0'	
2	'baselineEEG02.mat'	250	'Geodesic Sensor Net 64 2.0'	
3	'baselineEEG03.mat'	250	'Geodesic Sensor Net 64 2.0'	
4	'baselineEEG04.mat'	250	'HydroCel GSN 128 1.0'	
5	'baselineEEG05.mat'	250	'HydroCel GSN 128 1.0'	
6	'baselineEEG06.mat'	250	'HydroCel GSN 128 1.0'	
7	'baselineEEG07.mat'	250	'HydroCel GSN 128 1.0'	
8	'baselineEEG08.mat'	500	'HydroCel GSN 128 1.0'	
9	'baselineEEG09.mat'	500	'HydroCel GSN 128 1.0'	
10	'baselineEEG10.mat'	500	'HydroCel GSN 128 1.0'	
11				
12				

Workspace

Name ▲

Value

beapp_file_info_table

10x3 table

Information table for auditory (.mff) files:

Current Folder

Name ▾

beapp

set_beapp_path.m

set_beapp_def.m

prepare_to_run_main.m

beapp_main.m

beapp_gui.m

beapp_configure_settings.m

user_inputs

rerun_fselect_table.mat

beapp_userinputs.m

beapp_set_input_file_locations.m

beapp_file_info_table.mat

beapp_advinputs.m

run_templates

reference_data

Packages

functions

documentation

Variables – beapp_file_info_table

beapp_file_info_table

10x2 table

	1	2	3	4
	FileName	FileOffset		
1	'auditoryEEG01.mff'	0		
2	'auditoryEEG02.mff'	0		
3	'auditoryEEG03.mff'	0		
4	'auditoryEEG04.mff'	0		
5	'auditoryEEG05.mff'	0		
6	'auditoryEEG06.mff'	0		
7	'auditoryEEG07.mff'	8		
8	'auditoryEEG08.mff'	18		
9	'auditoryEEG09.mff'	18		
10	'auditoryEEG10.mff'	18		
11				
12				
13				

Workspace

Name ▲

Value

beapp_file_info_table

10x2 table